

Acetonitrile-d3

Art. ID TRC-A163892-5G
Unit 5 g
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

Description

Category: Intermediates, Miscellaneous Reagents /// Appearance: Clear colourless Oil /// Application Notes: Isotope labelled Acetonitrile is used as a polar aprotic solvent in organic synthesis. /// References: Hu, H.H., et al.: Phy. Chem. Phys., 15, 7666 (2013)| Kim, J.Y., et al.: J. Mater. Chem. A., et al.: 1, 5982 (2013)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Acetonitrile-d3	[2206-26-0]				