

Trace Metals - Ready-To-Use RM - Waste water (alternatively CRM see Code ERA-500)

Art. ID ERA-740
Unit 500 mL
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Aluminium (Al) | [7429-90-5] | 200 - 4000 | mg/L | | |
| | Antimony (Sb) | [7440-36-0] | 90 - 900 | mg/L | | |
| | Arsenic (As) | [7440-38-2] | 90 - 900 | mg/L | | |
| | Barium (Ba) | [7440-39-3] | 100 - 2500 | mg/L | | |
| | Beryllium (Be) | [7440-41-7] | 50 - 500 | mg/L | | |
| | Boron (B) | [7440-42-8] | 800 - 2000 | mg/L | | |
| | Cadmium (Cd) | [7440-43-9] | 100 - 1000 | mg/L | | |
| | Chromium (Cr) | [7440-47-3] | 100 - 1000 | mg/L | | |
| | Cobalt (Co) | [7440-48-4] | 100 - 1000 | mg/L | | |
| | Copper (Cu) | [7440-50-8] | 100 - 1000 | mg/L | | |
| | Iron (Fe) | [7439-89-6] | 200 - 4000 | mg/L | | |
| | Lead (Pb) | [7439-92-1] | 100 - 1500 | mg/L | | |
| | Manganese (Mn) | [7439-96-5] | 200 - 2000 | mg/L | | |
| | Molybdenum (Mo) | [7439-98-7] | 60 - 600 | mg/L | | |
| | Nickel (Ni) | [7440-02-0] | 200 - 2000 | mg/L | | |
| | Selenium (Se) | [7782-49-2] | 100 - 1000 | mg/L | | |
| | Silver (Ag) | [7440-22-4] | 100 - 1000 | mg/L | | |
| | Strontium (Sr) | [7440-24-6] | 50 - 500 | mg/L | | |
| | Thallium (Tl) | [7440-28-0] | 80 - 800 | mg/L | | |
| | Vanadium (V) | [7440-62-2] | 50 - 2000 | mg/L | | |
| | Zinc (Zn) | [7440-66-6] | 300 - 2000 | mg/L | | |