

N-Acetyl-DL-alanine

Art. ID TRC-A167965-5G
Unit 5 g
Deliverydetails No Dangerous Good

Description

Appearance: White to Off-White Solid /// Application Notes: N-Acetyl-DL-alanine is a derivative of DL-Alanine (A481400) and is used as a reagent in the synthesis of potent new agents for the treatment of epilepsy. ///
References: Conley, J.D., et al.: J. Med. Chem., 30, 567 (1987)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	N-Acetyl-DL-alanine	[1115-69-1]				