

N-Acetyl-DL-alanine-3,3,3-d3

Art. ID TRC-A167967-10MG
Unit 10 mg
Deliverydetails No Dangerous Good

Description

Application Notes: N-Acetyl-DL-alanine-3,3,3-d3 is an isotopic analog of N-Acetyl-DL-alanine (A167965), a derivative of DL-Alanine (A481400) and is used as a reagent in the synthesis of potent new agents for the treatment of epilepsy. /// References: Conley, J.D., et al.: J. Med. Chem., 30, 567 (1987)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	N-Acetyl-DL-alanine-3,3 ,3-d3					