

N-Acetyl-L-arginine Methyl Ester Hydrochloride

Art. ID TRC-A168360-100MG
Unit 100 mg
Deliverydetails No Dangerous Good

Description

Category: Amino Acids and Derivatives /// Application Notes: N-Acetyl-L-arginine Methyl Ester is used in the synthetic route to biologically important L-arginine peptides. /// References: Levison, P.R., et al.: Synthetic Communic., 13, 167 (1983)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	N-Acetyl-L-arginine Met hyl Ester Hydrochloride	[1784-05-0]				