

1-[(S)-3-(Acetylthio)isobutyryl]-L-proline Hydrate

Art. ID TRC-A188730-100MG
Unit 100 mg
Deliverydetails No Dangerous Good /not restricted

Description

Category: Miscellaneous Reagents /// Application Notes: 1-[(S)-3-(Acetylthio)isobutyryl]-L-proline Hydrate is used in the preparation of captopril derivatives, used as angiotensin converting enzymes. /// References: Li, H. et al.: Res. Chem. Intermed., 2, 621 (2013)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1-[(S)-3-(Acetylthio)isobutyryl]-L-proline Hydrate	[64838-55-7]				