

Adenosine 5'-Diphosphate Ammonium Salt Hydrate

Art. ID TRC-A281700-250MG
Unit 250 mg
Deliverydetails No Dangerous Good

Description

Category: Carbohydrates and Derivatives, Heterocycles, Nucleotides, Bases and Related Reagents, Pharmaceuticals, Intermediates and Fine Chemicals /// Appearance: Light Brown Solid /// Application Notes: A nucleoside diphosphate that interacts with a family of ADP receptors found on platelets (P2Y1, P2Y12 and P2X1), to induce platelet activation. ADP in the blood is converted to adenosine by the action of ecto-ADPases, inhibiting further platelet activation via adenosine receptors. Important in energy flow /// References: Gamba, G. et al.: Haematologia, 58, 1215 (1973)| Dejana, E. et al.: Thromb. Haemost., 39, 135 (1978)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Adenosine 5'-Diphosphat e Ammonium Salt Hydrate					