

## **Adenosine 5'-Triphosphate Disodium Salt Hydrate**

Art. ID TRC-A281810-50MG  
Unit 50 mg  
Deliverydetails Dry ice shipment /not restricted

### Description

Category: Carbohydrates and Derivatives, Heterocycles, Metabolites and Impurities, Nucleotides, Bases and Related Reagent, Phosphorylating and Phosphitylating Agents /// Appearance: White to Off-White Solid /// Application Notes: A metabolite of Adenosine (A280400), a multifunctional nucleoside triphosphate used in cells as a coenzyme of intracellular energy transfer. It transports chemical energy within cells for metabolism. /// References: Knowles, J., et al.: Annu. Rev. Biochem., 49, 877 (1980)

| Text/Information | Analyte/Parameter                               | CAS number   | Concentration/Value | Unit | Method | Source |
|------------------|-------------------------------------------------|--------------|---------------------|------|--------|--------|
|                  | Adenosine 5'-Triphosphate Disodium Salt Hydrate | [34369-07-8] |                     |      |        |        |