

Adipic Acid

Art. ID TRC-A291590-10G
Unit 10 g
Deliverydetails No Dangerous Good /not restricted

Description

Category: Miscellaneous Reagents /// Appearance: White Solid /// Application Notes: Adipic Acid is primarily used in the synthesis of nylon. It has been used as a reagent in the solid-state polymerization of nylon analogs. /// References: Kim, Y. et al.: Ind. Eng. Chem. Res., 51, 15801 (2012)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Adipic acid	[124-04-9]				