

Adipoyl Chloride

Art. ID TRC-A291700-5G
Unit 5 g
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN3265 II E2

Description

Category: Acyclics, Miscellaneous Reagents /// Appearance: Clear Colourless Oil /// Application Notes: Adipoyl Chloride can be used in the synthesis of nylon. Also it is used in the synthesis of chiral polymer for membrane application. /// References: Gaymans, R. et al.: J. Polym. Sci., 15, 535 (1977)| Hazarika, S. et al.: Polym. J., 41, 1067 (2009)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Adipoyl chloride	[111-50-2]				