

Alendronic Acid-d6 Sodium Salt Hydrate

Art. ID	TRC-A521257-2.5MG
Unit	2,5 mg
Deliverydetails	No Dangerous Good

Description

Category: Inhibitors, Isotope Labelled Compounds, Pharmaceuticals, Intermediates and Fine Chemicals, Phosphorylating and Phosphitylating Agents /// Appearance: White to Off-White Solid /// Application Notes: Alendronic Acid-d6 Sodium Salt is a labelled salt analogue of Alendronic Acid (A521250), a bisphosphonate used as a bone reabsorption inhibitor. /// References: Wang, X.M., et al.: Clin. J. Dent. Res., 3, 26 (2000)| Woo, T., et al.: Best Pract. Res. Clin. Rheumatol., 15, 469 (2001)| Benford, H.L., et al.: Bone, 28, 465 (2001)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Alendronic Acid-d6 Sodi um Salt Hydrate					