

D-allo-Threonine

Art. ID TRC-A548540-500MG
Unit 500 mg
Deliverydetails No Dangerous Good

Description

Category: Amines, Amino Acids and Derivatives, Chiral Reagents /// Application Notes: D-allo-Threonine is used to prepare serine and threonine beta-lactones as inhibitors of hepatitis A virus 3C cysteine proteinase. It is also used to synthesize the aglycon of ramoplanin A2 and ramoplanose with antibacterial activities. /// References: Lall, M., et al.: J. Org. Chem., 67, 1536 (2002)| Jiang, W., et al.: J. Am. Chem. Soc., 125, 1877 (2003)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	D-allo-Threonine	[24830-94-2]				