

Amidotrizoic Acid Dihydrate

Art. ID TRC-A576805-1G
Unit 1 g
Deliverydetails No Dangerous Good /not restricted

Description

Category: Aromatics, Pharmaceuticals, Intermediates and Fine Chemicals /// Appearance: White Solid ///
Application Notes: Amido-trizoic Acid Dihydrate is used as a density gradient reagent for blood cell separation.
Diagnostic aid (radiopaque medium). /// References: Kalmer, J.R., et al., J. Immunol. Methods, 110, 275 (1988),
Wang, Y., et al.: Br. J. Radiol., 170, 1229 (1997)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Amidotrizoic acid dihyd rate	[50978-11-5]				