

4-Amino-L-trans-proline

Art. ID TRC-A593885-250MG
Unit 250 mg

Description

Category: Amines, Chiral Reagents, Pharmaceuticals, Intermediates and Fine Chemicals /// Application Notes: 4-Amino-L-trans-proline is used to improve antagonistic activity of alpha-conotoxin analogues targeting alpha7 nicotinic acetylcholine receptors. /// References: Armishaw, C., et. al.: J. Biol. Chem., 284, 9498 (2009)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4-Amino-L-trans-proline	[16257-88-8]				