

DL-2-Aminocaprylic Acid

Art. ID	TRC-A603185-100MG
Unit	100 mg
Deliverydetails	No Dangerous Good /not restricted

Description

Category: Aliphatics /// Application Notes: DL-2-Aminocaprylic Acid is used as a reagent in the preparation of 1,5-disubstituted-2-aminoimidazoles which has antibiotic activity against a variety of bacterial strains and even multi drug resistant isolates. DL-2-Aminocaprylic Acid is also used as a reactant in the solution phase parallel synthesis of fenbufen (F247200) and ethacrynic acid (E676000). /// References: Harris, T.L., et al.: Bioorg. Med. Chem. Lett., 21, 4516 (2011)| Su, Y., et al.: Bioorg. Med. Chem. Lett., 21, 1320 (2011)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	DL-2-Aminocaprylic acid	[644-90-6]				