

(2R)-2-Amino-3-(4-aminophenyl)propanoic Acid

Art. ID TRC-A611345-100MG
Unit 100 mg
Deliverydetails No Dangerous Good

Description

Application Notes: (2R)-2-Amino-3-(4-aminophenyl)propanoic Acid is an amino acid that is used as a reactant in the development of human calcitonin gene-related peptide (CGRP) receptor antagonists for treatment of acute migraine. /// References: Rudolf, K., et al.: J. Med. Chem., 48, 5921 (2005)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(2R)-2-Amino-3-(4-amino phenyl)propanoic Acid	[102281-45-8]				