

p-Amino-L-phenylalanine

Art. ID TRC-A622376-2.5G
Unit 2,5 g
Deliverydetails No Dangerous Good

Description

Category: Amino Acids and Derivatives, Aromatics, Chiral Reagents, Pharmaceuticals, Intermediates and Fine Chemicals /// Appearance: Beige Solid /// Application Notes: p-Amino-L-phenylalanine is an analog of L-Phenylalanine (P319415). /// References: Emes, A., et al.: Can. J. Microbiol., 20, 347 (1974), Green, J., et al.: J. Biol. Chem., 266, 12971 (1991), Chang, Z., et al.: Microbiology, 147, 2113 (2001)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	p-Amino-L-phenylalanine	[943-80-6]				