

4-Aminophenylmercuric Acetate

Art. ID TRC-A626025-250MG
Unit 250 mg
Deliverydetails De Minimis UN2025 II E4

Description

Appearance: Solid /// Application Notes: 4-Aminophenylmercuric Acetate used in combination with heparin is said to be effective in activating matrix metalloproteinase-2 (MMP-2) in peripheral blood of healthy humans, according to biological studies. /// References: Mannello, F., et al.: Clin. Biochem., 41, 1466 (2008)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4-Aminophenylmercuric a cetate	[6283-24-5]				