

4-Aminohippuric-d4 Acid

Art. ID TRC-A627972-2.5MG
Unit 2,5 mg
Deliverydetails No Dangerous Good

Description

Category: Amines, Aromatics, Diagnostics, Isotope Labelled Compounds /// Appearance: Off-White to Pale Yellow Solid /// Application Notes: Labelled 4-Aminohippuric Acid is used in the measurement of renal plasma flow as a diagnostic tool which may be applied towards kidney disorders. /// References: Phillips, R. et al.: Am. J. Physiol., 152, 523 (1948)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4-Aminohippuric-d4 Acid	[1219805-41-0]				