

O-Acetyl-L-homoserine Hydrochloride, 95%

Art. ID TRC-A178300-10MG
Unit 10 mg

Description

Category: Amino Acids and Derivatives /// Appearance: Off-White to Pale Yellow Solid /// Application Notes: A versatile synthon for the synthesis of L-Homoserine peptides and 3-Amino-2-pyrrolidinones. /// References: Baier, M., et al.: Plant Physiol., 145, 1600 (2007), XXX

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	O-Acetyl-L-homoserine H ydrochloride, 95%	[250736-84-6]				