

2-Aminopropane-1,3-diol

Art. ID TRC-A628115-1G
Unit 1 g
Deliverydetails De Minimis UN3263 II E2

Description

Category: Amines, Miscellaneous Reagents /// Appearance: White to Off-White Solid /// Application Notes: 2-Aminopropane-1,3-diol is used in the synthesis of novel potent inhibitors of human uridine phosphorylase used to reduce toxicity of 5-fluorouracil without affecting its antitumor activity. /// References: Renck, D. et al.: J. Med. Chem., 56, 8892 (2013)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Serinol	[534-03-2]				