

N-Acetyl-3-iodo-L-tyrosine Amide

Art. ID TRC-A179410-250MG
Unit 250 mg

Description

Category: Amino Acids and Derivatives, Aromatics, Chiral Reagents /// Appearance: Off-White Solid /// Application Notes: Used in the synthesis of 3-iodo-L-thyronine and its iodinated derivatives. /// References: Block, P., et al.: J. Med. Chem., 19, 1067 (1976), XXX

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	N-Acetyl-3-iodo-L-tyrosine Amide	[59302-19-1]				