

## **DL-Alaninol**

Art. ID TRC-A628760-5G  
Unit 5 g  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN2735 II E2

### Description

Appearance: Clear Colourless Liquid /// Application Notes: DL-Alaninol is a chemical reagent used in organic synthesis. Used in the preparation of polyunsaturated fatty acid amides showing anti-proliferative properties. /// References: Tremblay, H. et al.: Bioorg. Med. Chem. Lett., 24, 5635 (2014)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	DL-Alaninol	[6168-72-5]				