

## **2-Aminopyrimidine-5-boronic Acid Pinacol Ester**

Art. ID TRC-A629055-500MG  
Unit 500 mg  
Deliverydetails No Dangerous Good

### Description

Category: Aromatics, Boron Derivatives, Heterocycles, Intermediates, Miscellaneous Reagents /// Appearance: Light Beige Solid /// Application Notes: 2-Aminopyrimidine-5-boronic Acid Pinacol Ester is used in the preparation of phosphoinositide-3-kinases (PI3K) inhibitors. /// References: Burger, M.T. et al.: ACS Med. Chem. Lett., 2, 774 (2011)| Heffron, T.P. et al.: J. Med. Chem., 55, 8007 (2012)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2-Aminopyrimidine-5-boronic Acid Pinacol Ester	[402960-38-7]				