

rac-2-Aminobutyric Acid-d3

Art. ID TRC-A602937-25MG
Unit 25 mg

Description

Category: Amino Acids and Derivatives, Chiral Reagents, Isotope Labelled Compounds /// Appearance: White to Off-White Solid /// Application Notes: A descriptor amino acid QSARs peptide. /// References: Kidera, A., et al.: J. Protein Chem., 4, 23

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	rac-2-Aminobutyric Acid -d3	[1219373-19-9]				