

## **L-Arabinitol**

Art. ID TRC-A735610-2.5G  
Unit 2,5 g  
Deliverydetails No Dangerous Good /not restricted

### Description

Category: Carbohydrates and Derivatives, Diagnostics /// Appearance: Off-White Solid /// Application Notes: L-Arabinitol is the stereoisomer of D-Arabinitol (A735600), both used as diagnostic tools for the detection of fungal infections in neutropenic patients. /// References: Salonen, J.H., et al.: Eur. J. Clin. Microbiol. Infect. Dis., 20, 179 (2001)| Sigmundsdottir, G. et al.: J. Clin. Microbiol., 38, 3039 (2000)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	L(-)-Arabit	[7643-75-6]				