

Arachidic Acid

Art. ID TRC-A764990-100MG
Unit 100 mg
Deliverydetails No Dangerous Good /not restricted

Description

Category: Aliphatics, Glycerols, Fatty Acid Derivatives and Lipids /// Appearance: White to Off-White Solid ///
Application Notes: Arachidic acid is a saturated fatty acid and a minor constituent of peanut oil, coconut oil and corn oil. Arachidic acid is used for the production of detergents, photographic materials and lubricants. /// References: Kumar, S. R., et al.: J Chemtracks., 13, 277 (2011)| Poyarkova, T.N., et al: J. App. Chem. 85, 1259 (2012)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arachidic acid, No. of	[506-30-9]				
	C: C20:0					