

Aripiprazole Dimer

Art. ID TRC-A771025-2.5MG

Unit 2.5 mg

Description

Category: Aromatics, Heterocycles, Pharmaceuticals, Intermediates and Fine Chemicals /// Appearance: Tan to Light Beige Solid /// Application Notes: Aripiprazole Dimer is an impurity of Aripiprazole (A771000), a selective dopamine D2-receptor antagonist with dopamine autoreceptor agonist activity. /// References: Canive, J.M., et al.: Psychopharmacol. Bull., 34, 101 (1998)| Oshiro, Y., et al.: J. Med. Chem., 41, 658 (1998)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aripiprazole Dimer	[1797986-18-5]				