

A77 1726 (E/Z) Mixture

Art. ID TRC-A771720-50MG
Unit 50 mg
Deliverydetails No Dangerous Good /not restricted

Description

Category: Metabolites and Impurities, Tyrosine Kinase Inhibitors, Pharmaceuticals, Intermediates and Fine Chemicals /// Appearance: White Solid /// Application Notes: rac A77 1726 is the active metabolite of Leflunomide, a potent disease-modifying antirheumatic drug used in the treatment of rheumatoid arthritis. LEF-M interferes with dendritic cell function. The second Phase III study evaluating the efficacy of Genzyme's teriflunomide therapy against relapsing forms of multiple sclerosis found the once-daily oral drug to be no less effective than injectable interferon beta 1a (Rebif®). /// References: Kirsch, B.M., et al.: Arthritis Research and Therapy, 7, 3, R694, (2005)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2-Cyano-3-hydroxy-N-(4'-trifluoromethylphenyl)-crotonamide	[108605-62-5]				