

## **Arsenic Oxide (As<sub>2</sub>O<sub>3</sub>)**

Art. ID TRC-A774130-10G  
Unit 10 g  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1561 II E4

### **Description**

Category: Miscellaneous Reagents /// Appearance: White to Pale Yellow Solid /// Application Notes: Arsenic Oxide ((As<sub>2</sub>O<sub>3</sub>)) is used as a precursor to forestry products, in colorless glass production, and in electronics. In combination with copper(II), acetate arsenic trioxide gives Copper Acetoarsenite (C685255) as a rodenticide. /// References: Grund, S., et al.: Ullmann's Encyclopedia of Industrial Chemistry, Weinheim: Wiley-VCH (2005)| Gaines, T.: Toxicol. Appl. Pharm., 2, 88 (1960)| Jameson, F.: Paint Manu., 13, 6 (1943)

| Text/Information | Analyte/Parameter   | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|---------------------|-------------|---------------------|------|--------|--------|
|                  | Arsenic (III) oxide | [1327-53-3] |                     |      |        |        |