

AS-252424

Art. ID TRC-A780000-5MG
Unit 5 mg
Deliverydetails No Dangerous Good /not restricted

Description

Category: Amines, Aromatics, Heterocycles, Inhibitors, Sulfur and Selenium Compounds /// Appearance: Yellow Solid /// Application Notes: AS-252424 is a novel, potent, and selective phosphoinositide 3-kinase gamma (PI3Kgamma) inhibitor. It is also an inhibitor of Casein kinase and PI3K p110alpha. /// References: Pomel, V., et al.: J. Med. Chem., 49, 3857 (2006)| Edling, C.E., et al.: Clin. Cancer Res., 16, 4928 (2010)| Condliffe, A.M., et al.: Blood, 106, 1432 (2005)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	AS-252424	[900515-16-4]				