

trans-alpha-Asarone

Art. ID TRC-A780205-500MG
Unit 500 mg
Deliverydetails No Dangerous Good /not restricted

Description

Category: Aromatics, Miscellaneous Reagents /// Appearance: White to Off-White Solid /// Application Notes: trans-alpha-Asarone is used in the synthesis of anti-inflammatory neolignan. Also used in lignan syntheses for antiallergic agents. /// References: Kadota, S. et al.: Tetrahedron Lett., 28, 2857 (1987)| Mori, K. et al.: Tetrahedron, 42, 523 (1986)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	alpha-Asarone	[2883-98-9]				