

Monomethyl Arsonothioic Acid (80%)

Art. ID	TRC-A788260-1MG
Unit	1 mg
Deliverydetails	No Dangerous Good

Description

Appearance: Off-White Solid /// Application Notes: As-Methyl Arsonothioic Acid is used as a standard for arsenic inspection. It is used as an aid to help in determining the level of toxicity in humans. Metabolite in some mammals and microbiota. /// References: Naranmandura, H. et al.: Chem. Res. Toxicol., 24, 1586 (2011)| Bu, N. et al.: Metallom., 3, 1064 (2011)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Monomethyl Arsonothioic Acid (80%)	[937022-98-5]				