

L-Aspartic Anhydride Hydrochloride

Art. ID TRC-A788780-500MG
Unit 500 mg
Deliverydetails No Dangerous Good

Description

Application Notes: L-Aspartic Anhydride Hydrochloride is a chemical acylating reagent used in friedel-crafts reactions such as in the preparation of L-homophenylalanine. /// References: Lin, W. et al.: Synthesis, 7, 1007 (2001)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	L-Aspartic Anhydride Hydrochloride	[34029-31-7]				