

L-2-Amino-3-[thiomethyl-1-(1-oxyl-2,2,5,5-tetramethyl-3-pyrrolin-3-yl)]propanoic Acid Dihydrochloride

Art. ID TRC-A630400-5MG

Unit 5 mg

Description

Category: Amino Acids and Derivatives, Spin Labeling Compounds, Sulfur and Selenium Compounds ///

Appearance: Light Yellow to Dark Orange Solid /// Application Notes: A synthetic spin labelled amino acid that can be used as a biophysical probe to detect structural changes in proteins and the interaction of proteins with other macromolecules. /// References: Cornish, V.W., et al.: Proc. Natl. Acad. Sci. USA, 91, 2910 (1994)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	L-2-Amino-3-[thiomethyl -1-(1-oxyl-2,2,5,5-tetr amethyl-3-pyrrolin-3-y l)]propanoic Acid Dihyd rochloride					