

Androsterone-D5

Art. ID TRC-A637537-1MG
Unit 1 mg

Description

Category: Isotope Labelled Compounds, Metabolites and Impurities, Pharmaceuticals, Intermediates and Fine Chemicals, Steroids, /// Application Notes: A labelled steroid hormone with weak androgenic activity. It is metabolized from Testosterone (T155000) in the liver and is used in doping analysis to detect testosterone misuse. Controlled Substance. /// References: Park, J., et al.: J. Anal. Toxicol., 14, 66 (1990),| Palonek, E., et al.: J. Steroid Mol. Biol., 55, 121

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Androsterone-D5					