

5-Azido Uridine

Art. ID TRC-A926800-25MG
Unit 25 mg

Description

Category: Nucleotides, Bases and Related Reagents, Pharmaceuticals, Intermediates and Fine Chemicals ///
Appearance: Yellow Solid /// Application Notes: Used in a method and kit for labeling nucleic acid in living cells. The labeled nucleic acid can be used for detecting cell proliferation and cell apoptosis, labeling plasmid vector, virus, small interfering RNA and microRNA, and studying early development.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|----------------|---------------------|------|--------|--------|
| | 5-Azido Uridine | [1355028-82-8] | | | | |