

Aripiprazole-d8

Art. ID TRC-A771002-1MG
Unit 1 mg

Description

Category: Heterocycles, Isotope Labelled Compounds, Neurochemicals, Pharmaceuticals, Intermediates and Fine Chemicals /// Appearance: Off-White to Pale Beige Solid /// Application Notes: Aripiprazole-d8 is a deuterated version of Aripiprazole, a selective dopamine D2-receptor antagonist with dopamine autoreceptor agonist activity. /// References: Canive, J.M., et al.: Psychopharmacol. Bull., 34, 101 (1998), Oshiro, Y., et al.: J. Med. Chem., 41, 658 (1998)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aripiprazole-d8	[1089115-06-9]				