

## **Aripiprazole-d8 (Butyl-d8)**

Art. ID TRC-A771003-1MG  
Unit 1 mg  
Deliverydetails De Minimis

### Description

Category: Heterocycles, Isotope Labelled Compounds, Neurochemicals, Pharmaceuticals, Intermediates and Fine Chemicals /// Appearance: White Solid /// Application Notes: Aripiprazole-d8 (Butyl-d8) is deuterated version of Aripiprazole, a selective dopamine D2-receptor antagonist with dopamine autoreceptor agonist activity. /// References: Canive, J.M., et al.: Psychopharmacol. Bull., 34, 101 (1998), Oshiro, Y., et al.: J. Med. Chem., 41, 658 (1998)

| Text/Information | Analyte/Parameter          | CAS number     | Concentration/Value | Unit | Method | Source |
|------------------|----------------------------|----------------|---------------------|------|--------|--------|
|                  | Aripiprazole-d8 (Butyl-d8) | [1089115-04-7] |                     |      |        |        |