

D-Azetidine-2-carboxylic Acid

Art. ID TRC-A812480-500MG
Unit 500 mg

Description

Category: Amino Acids and Derivatives /// Appearance: White to Off-White Solid /// Application Notes: A four-membered ring analog of L-Proline. A useful intermediate in the synthesis of polypeptides. /// References: Barber, M., et al.: Int. J. Pept. Protein Res., 14(3), 247 (1979), Tsai., F.H., et al.: Biopolymers, 30(11-12), 1039 (1990)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	D-Azetidine-2-carboxyli c Acid	[7729-30-8]				