

Heavy Metals in Soft Drink Reference Material (Material: TET028RM)

Art. ID FAPAS-FCCM3-DRN29RM
Unit approx. 50 mL
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Antimony (Sb)	[7440-36-0]	239,4 ± 7,5	µg/L		
	Arsenic (As)(total)		907,5 ± 42,6	µg/L		
	Cadmium (Cd)	[7440-43-9]	246,0 ± 7,2	µg/L		
	Chromium (Cr)	[7440-47-3]	751,1 ± 42,4	µg/L		
	Copper (Cu)	[7440-50-8]	415,9 ± 26,4	µg/L		
	Zinc (Zn)	[7440-66-6]	940,5 ± 62	µg/L		