

o-Acetyltoluene

Art. ID TRC-A189000-25G
Unit 25 g
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ)

Description

Category: Aromatics, Miscellaneous Reagents /// Appearance: Clear Colourless Oil /// Application Notes:
o-Acetyltoluene can be used as a reagent to prepare diaryl dihydropyrazoles as inhibitors of mitotic kinesin spindle protein (1). It is also used to synthesize aryl- and indolyl-cyanoaminopyridines with anti-inflammatory, analgesic, and antipyretic activities (2).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	o-Acetyltoluene	[577-16-2]				