

Adipic Acid Dihydrazide

Art. ID TRC-A291595-50G
Unit 50 g
Deliverydetails No Dangerous Good /not restricted

Description

Category: Amines, Miscellaneous Reagents /// Appearance: White Solid /// Application Notes: Adipic Acid Dihydrazide is used to functionalize magnetic nanoparticles for glycopeptide enrichment and identification. Adipic Acid Dihydrazide can also be used as a chain extension for liquid rubber. /// References: Cao, Q., et. al.: Analyst, 139, 603 (2014)| Cutts, E.: Chem. Industry, 7, 280 (1974)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Adipic acid dihydrazide	[1071-93-8]				