

## **Adipoyl Chloride**

Art. ID TRC-A291700-25G  
Unit 25 g  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN3265 II E2

### **Description**

Category: Acyclics, Miscellaneous Reagents /// Appearance: Clear Colourless Oil /// Application Notes: Adipoyl Chloride can be used in the synthesis of nylon. Also it is used in the synthesis of chiral polymer for membrane application. /// References: Gaymans, R. et al.: J. Polym. Sci., 15, 535 (1977)| Hazarika, S. et al.: Polym. J., 41, 1067 (2009)|

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|------|--------|--------|
|                  | Adipoyl chloride  | [111-50-2] |                     |      |        |        |