

L-Alanyl-L-proline

Art. ID TRC-A511340-500MG
Unit 500 mg
Deliverydetails /not restricted

Description

Category: Amines, Amino Acids and Derivatives, Heterocycles, Pharmaceuticals, Intermediates and Fine Chemicals
/// Appearance: White Solid /// Application Notes: L-Alanyl-L-proline is a dipeptide that contains a sterically constrained proline, a well-known turn inducer in proteins. /// References: Moriuchi, T., et al.: J. Organomet. Chem., 589, 50 (1999)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Ala-Pro hydrate	[13485-59-1]				