

Allyltrimethylsilane

Art. ID TRC-A572350-10G
Unit 10 g
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1993 II E2

Description

Category: Miscellaneous Reagents /// Appearance: Clear Colourless Oil /// Application Notes: Allyltrimethylsilane is used in the allylation of aldehydes, imines, allylic and benzylic alcohols, and chiral alpha-keto-amides that are derived from (S)-proline esters. /// References: Yadav, J.S., et al.: Tetrahedron Lett., 43, 3783 (2002)| Aggarwal, V.K., et al.: Tetrahedron Lett., 37, 3745 (1996)| Wang, D.K., et al.: J. Org. Chem., 64, 4233 (1999)| Kaur, G., et al.: Tetrahedron Lett., 38, 2521 (1997)| Soai, K., et al.: J. Chem. Soc., Chem. Commun., 15, 1016 (1984)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Allyltrimethylsilane	[762-72-1]				