

4-Amino-3,5-bis(4-pyridyl)-1,2,4-triazole

Art. ID TRC-A580850-100MG

Unit 100 mg

Description

Application Notes: 4-Amino-3,5-bis(4-pyridyl)-1,2,4-triazole can be used to prepare cadmium aminoisophthalate pyridinyltriazolamine metal-organic frameworks for detection applications. It has also been utilized to synthesize a linear trinuclear Cu(II) complex with magnetic properties. /// References: Wang, H., et al.: RSC Adv. 5, 48881 (2015)| Zhang, Q., et al.: J. Coord. Chem. 64, 2110 (2011)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4-Amino-3,5-bis(4-pyridyl)-1,2,4-triazole	[38634-05-8]				