

2-Amino-6-methylpyridine

Art. ID TRC-A613985-25G
Unit 25 g
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN2811 III E1

Description

Application Notes: 2-Amino-6-methylpyridine is a reagent used in the identification of selective inhibitor for phosphatidylinositol 3-kinase (PI3Kalpha). /// References: Miller, M. S., et al.: Bioorg. Med. Chem., 25, 1481-1486 (2017)|

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2-Amino-6-methylpyridin	[1824-81-3]				
	e					