

DL-2-Amino-4-phosphonobutanoic acid (AP4)

Art. ID TRC-A627558-100MG
Unit 100 mg
Deliverydetails No Dangerous Good

Description

Application Notes: DL-2-Amino-4-phosphonobutanoic acid (AP4) is an amino-acid derivative reported to block cone signals in the rat retina. Improves learning and memory processes in passive and active avoidance situations in rats. /// References: Zalewska, A. et al.: Polish J. Pharmacol., 47, 9 (1995)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	DL-2-Amino-4-phosphonob	[6323-99-5]				
	utanoic acid (AP4)					