

4-Aminohippuric Acid

Art. ID TRC-A627970-100G
Unit 100 g
Deliverydetails No Dangerous Good /not restricted

Description

Category: Amines, Aromatics, Diagnostics /// Appearance: White Solid /// Application Notes: 4-Aminohippuric Acid is used in the measurement of renal plasma flow as a diagnostic tool which may be applied towards kidney disorders. /// References: Phillips, R. et al.: Am. J. Physiol., 152, 523 (1948)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aminohippuric Acid	[61-78-9]				